



## New Image Shape-to-Fit FormaFlex Skin Barrier with Tape 1410X

**Series:** 1410X



### Product Benefits

- A Shape-to-Fit Skin Barrier can be stretched to the shape and size of the stoma. This skin barrier does not require the use of scissors. A possible choice when the stoma is not round or still changing sizes.
- Confidence assured with the Floating Flange; provides utmost security with the patented locking system that audibly "clicks" into place.
- Eliminates pressure on tender post operative abdomens with the exclusive Floating Flange.

### Product Features

- FormaFlex Skin Barrier
- Color Matching System
- Floating Flange
- Tape Border
- Shape-to-fit

### Products in this Series

Stock #	Barrier Opening	Color ID	Flange Size	Unit of Measure
14102	Up to 1-1/4" (32mm)	Green	1-3/4" (44mm) Flange	Box of 5
14103	Up to 1-11/16" (43mm)	Red	2-1/4" (57mm) Flange	Box of 5
14104	Up to 2-1/4" (57mm)	Blue	2-3/4" (70mm) Flange	Box of 5

**HCPCS Code:** A4409

